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Got rbST in your milk?

Dairy co-op bows to pressure to stop use of hormone

George Raine, Chronicle Staff Writer
Sunday, March 25, 2007

In January, California Dairies' board voted to ask its members not to inject synthetic bovine growth hormone into their cows. If they do, their milk will have to be segregated and they'll pay a surcharge.


Consumer demand is obvious," Cotta said.

Starbucks spokeswoman Sanja Gould in Seattle said the company is already rbST-free in stores in Northern California, New England, New Mexico, Montana, Oregon, Washington, Idaho and Alaska. That represents 37 percent of its dairy volume. The company will be free of bovine growth hormone in other markets as supply becomes available, Gould said.

Recombinant bovine somatotropin, or rbST, was approved by the Food and Drug Administration 14 years ago. Injected every two weeks into cows, it sustains lactation by stimulating cows' appetites so they eat more and produce more milk, perhaps an extra 5 quarts per day.

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


exist regarding its use, but the lack of confirmatory methods [1] for its detection makes it difficult to apply these regulations. It turns out to be an international issue in terms of animal doping and in food safety as well. Indeed, residues of rbST can be present in the milk of dairy animals treated with this hormone.

In order to detect residues of rbST in milk, the choice has been made to focus the analysis on the tryptic N-terminal peptide of the protein, specifically the difference between the endogenous and recombinant forms. The N-terminal amino acid alanine that is present in the endogenous form is replaced by a methionine in the recombinant one [2].

Abstract

Recombinant bovine somatotropin (rbST), also called growth hormone, is a protein hormone used in dairy farming to enhance milk production. A method has been developed for the detection of rbST in milk by ESI(+)-LC-MS/MS. This method allowed a detection limit of 20 pg of tryptic N-terminal peptide rbST in standard solution injected on-column and was successfully applied to extracts obtained from milk samples spiked with 50 ng/mL (2.3 pmol/mL) rbST.

This application describes a method for the detection of rbST by ESI(+)-LC-MS/MS. The method was successfully applied to extracts from milk samples spiked with rbST.

Experimental

Standards of Proteins and Peptides

Protein standards of rbST and recombinant equine somatotropin, reST (EquiGen-5), were obtained

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
Simultaneous Quantitative Screening and Qualitative Confirmation of 300 Pesticides Using the New 3200 Q TRAP® LC/MS/MS System

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
Automated On-line Solid-Phase Extraction–Liquid Chromatography–Electrospray Tandem Mass Spectrometry Method for the Determination of Ochratoxin A in Wine and Beer

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Department of Chemistry, "La Sapienza" University, Piazzale Aldo Moro 5, 00185 Roma, Italy

Ochratoxine A
(PM = 403)

Liquid chromatography/tandem mass spectrometric confirmatory method for determining aflatoxin M1 in cow milk
Comparison between electrospray and atmospheric pressure photoionization sources

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
Department of Chemistry, "La Sapienza" University, Piazzale Aldo Moro 5, 00185 Roma, Italy
Received 13 July 2005; received in revised form 11 September 2005; accepted 25 September 2005

Masked Mycotoxins: Determination of a Deoxynivalenol Glucoside in Artificially and Naturally Contaminated Wheat by Liquid Chromatography–Tandem Mass Spectrometry

FRANZ BERTHILLER,¹ CHIARA DALL'ASTA,¹ RAINER SCHUEHMACHER,^{2*} MARC LEMMENS,² GERHARD ADAM,² AND RUDOLF KRSKA¹

Christian Doppler Laboratory for Mycotoxin Research and Institute for Plant Production Biotechnology, Department for Agrobiotechnology (IFA-Tulln), University of Natural Resources and Applied Life Sciences, Vienna, Konrad Lorenz Strasse 20, A-3430 Tulln, Austria, and Institute of Applied Genetics and Cell Biology, Department of Applied Plant Sciences and Plant Biotechnology, University of Natural Resources and Applied Life Sciences, Vienna, Muthgasse 18, A-1190 Vienna, Austria

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



5ppt (pg/mL) on API 3200™ LC/MS/MS system

Malachite Green Testing Summary

- Detection of Malachite green and Leucomalachite green by LC/MS/MS in less than 10min
- API 3200™ LC/MS/MS system's data
 - LOD Malachite green (quantifier MRM) < 5ppt
 - Linear range Malachite green (both MRMs): 5-10000ppt
 - LOD Leucomalachite green (quantifier MRM) < 5ppt
 - Linear range Leucomalachite green (both MRMs): 50-100000ppt

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Characterization of acylated flavonoid-O-glycosides and methoxylated flavonoids from *Tagetes maxima* by liquid chromatography coupled to electrospray ionization tandem mass spectrometry

Irene Parejo *et al.*

Departament de Productes Naturals, Universitat de Barcelona, Spain

Quantitative determination of Astragaloside IV, a natural product with cardioprotective activity, in plasma, urine and other biological samples by HPLC coupled with tandem mass spectrometry

Weidong Zhang^{a,b,*}, Chuan Zhang^a, Runhui Liu^a, Huiliang Li^a,
Jintao Zhang^c, Chen Mao^c, Chunlin Chen^c

Simultaneous determination of four active alkaloids from a traditional chinese medicine *Corydalis saxicola* Bunting. (Yanhuanglian) in plasma and urine samples by LC-MS-MS

Hui-Liang Li^a, Wei-Dong Zhang^{a,b,*}, Run-Hui Liu^a, Chuan Zhang^a, Ting Han^d,
Xiang-Wei Wang^c, Xiao-Lin Wang^c, Jian-Bao Zhu^c, Chun-Lin Chen^c

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